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Saphire Bond

Light cure single component adhesive

Saphire Bond is a simple to use **one bottle light cure single component priming and dentinoenamel adhesive in one**. It is designed for strong bonding of **composites, compomers** and metals to **enamel and dentin and nonprecious and precious metals**. Saphire Bond can be also used for priming the root canals before filling of the root canals or cementation of endodontic posts.

Saphire Bond is further designed as a primer for indirect adhesive restorations, e.g. ceramic, metal and composite inlays, onlays, veneers, crowns and bridges with a self or dual curing composite cement.

Its strong adhesion to etched enamel or dentin proceeds on principles similar to that occurring with glassionomer cements. Good, long lasting adhesive strength and good biocompatibility are attained by methacrylate grafted polycarboxylic acids.

Saphire Bond is compatible with all current brands of visible light cure composite restorative materials. It is ethanol-based and hydrophilic. Therefore Saphire Bond can be used on slightly moist dentin surfaces (wet bonding technique).

In the rare cases a dual or self cure priming and bonding system is recommended Saphire Bond can get dual cure which is applied as a 1:1-mixture of Saphire Bond and Saphire Bond Activator. The activator is available separately. For these applications refer to recommendations for use of Saphire Bond Activator.

1. Bonding of Light Cure Composite and Composer Filling Restorations

1.1. Isolation

Rubber dam is the recommended method of isolation.

1.2. Cavity Preparation

Clean the tooth with flour of pumice and water prior to preparation. Prepare the cavity with minimal tooth reduction. Margins should have a slight (**0.5 - 1.0 mm**) bevel placed in the enamel to increase the surface area for greater bond strength.

1.3. Pulp Protection

Cavity floor of deep excavations should be covered with a thin layer of hard setting calciumhydroxide material.

1.4. Enamel and Dentine Conditioning

Recommended is the total etch technique. Apply **Saphire Conditioner** onto the enamel and dentine surfaces beginning with the enamel bevels. Leave the etching gel in place for **20 seconds**. On primary teeth and teeth high in fluoride, a **60 second** etch is recommended. Rinse for 20 seconds with water. Dry it in a water and oil-free airstream, but **do not** desiccate. A slightly wet dentine surface is important for the function of Saphire Bond. The etched enamel bevel should have a chalky white appearance.

Etching precaution:

It is essential, that etched areas are not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water.

For the typical indications of compomers etching is mostly not recommended, but it improves considerably the bonding with tooth material.

1.5. Application of Saphire Bond

Apply Saphire Bond generously with a brush onto the enamel and dentin surfaces for **30 seconds** with agitation. Dry cautiously with oil free air for 15 seconds. Do not desiccate the dentin.

Cure the Saphire Bond coating by exposing its entire area to a dental halogen light unit for **20 seconds** before application of a second layer of Saphire Bond.

Notes: Do not rinse off the Saphire Bond! If not used immediately, place dispensed Saphire Bond in subdued light to prevent premature polymerization by incident light. The Saphire Bond will **not** selfcure.

1.6 Application of a Second Layer of Saphire Bond

Apply again the Saphire Bond generously with a brush onto the adhesive surfaces for **30 seconds** with agitation.

Thoroughly dry again the primed areas with oil free air for **15 seconds** and cure the Saphire Bond coating by exposing its entire area to a dental halogen light unit for **20 seconds** before placement of a light cure composite.

1.7. Restorative Placement, Cure and Finishing

Refer to manufacturers instructions for placement, curing and finishing of **light cure restorative materials**.

2. Sealing of the Cavity before Amalgam Filling Restorations

2.1. Application of Light Cure Saphire Bond to Dentin/Enamel

After isolation of tooth, preparation of cavity, and in case of need pulp protection according point 1.1., 1.2. and 1.3., apply Saphire Bond **in two layers** by the same way as in **1.5** and **1.6** described.

It is essential that the primed areas are dry and contaminant free for the amalgam application. **Because light curing under amalgam restorations is impossible, and the primer does not selfcure** it is essential to cure the coating by exposing its entire area to a dental halogen light unit for **20 seconds before** application of the amalgam filling. Alternatively, please use the Saphire Bond in the dual cure mode (in addition with Saphire Bond Activator).

2.2. Application and Condensing of Amalgam

Apply amalgam according instructions of the manufacturer.

3. Priming of root canals

3.1. Preparation

Use of a rubber dam to isolate the tooth is strongly recommended.

Prepare and clean the root canal with hydrogen peroxide (3%) or sodium hypochlorite solution.

For fixing of endodontic posts with composite cements: Following the directions of the selected post manufacturer, prepare the canal to accept the appropriate length post. Trial fit the post and adjust for length, if necessary. Remove the post and prepare the canal for the endodontic post cementation.

3.2. Priming of root canals

Because the surfaces of the root canal surfaces are not accessible by a brush, use a syringe for the application of Saphire Bond and flush it into the prepared root canals. The flushing must go on over a time of **60 seconds** discontinuous in **4 cycles**. All the dentin surfaces must kept wet with the primer over the specified time by generously brushing the accessible surfaces with a brush.

An excess of the primer in the canals is to be removed with soft tissue paper. Dry the root canal cautiously with oilfree air for 15 seconds and cure the Saphire Bond coating with a dental halogen light unit for **20** seconds.

3.3. Cement placement

Refer to manufacturer instructions for placement and curing of the composite cements.

4. Application of Saphire Bondon Precious and Nonprecious Metals and Ceramics for repairing of crowns and bridges

4.1. Treatment of the surfaces of the prosthetic works

Rough surfaces of metallic or fractured ceramic surfaces with a diamond bur (40µm) or with suitable intraoral sand blasting. Treat the roughen surfaces with Saphire Conditioner, then rinse with plenty of water and dry afterwards with oil and waterfree air. Prevent contact of the etchant with skin eyes and gingival tissue.

Optional for ceramics:

If the prosthetic work can be done extraorally, ceramics can be etched with 1% hydrofluoric acid, but caution: hydrofluoric acid is strong corrosive. Do never use it intraorally.

It is recommended to prime the prosthetic work first with a suitable silane primer prior to application of Saphire Bond. Refer to instruction for use of the silane primer.

4.2. Application of Saphire Bond

Apply Saphire Bond generously with a brush onto the metallic surfaces and dry these surfaces with air that is oil and water-free for 15 seconds. Then cure the primer with a dental halogen light unit for **20 seconds** and start the repair work.

5. Reparations of fillings with amalgam

5.1. Treatment of the amalgam surfaces

Rough metallic surfaces with a diamond bur (40µm) or with suitable intraoral sand blasting.

Areas of enamel & dentine are to etch by the same way as in **1.4.** described

5.2. Application of Saphire Bond

Apply Saphire Bond in two layers according Chapter **1.5.** and **1.6.**

5.3. Application of a composite

For repairing the amalgam filling, applicate a light curing composite according manufacturers instructions.

6. Storage of Saphire Bond

Do not store above 25 °C (78 °F)! Avoid storage in direct sunlight. Do not use after expiration date.

7. Technical Data

Time for light cure with a dental halogen unit	20 sec
Setting-time of the 1:1-mix of Saphire Bond / Saphire Bond Activator for autocure (without any light; air excluded)	ca. 3 min

(refer to recommendation for use for Saphire Bond Activator).